

# LOCAL FIELD POTENTIAL -LFP RECORDING




# LFP CAPABILITIES

What is LFP

Measure the overall activity of an area rather than individual neurons

Differentiate borders between quiet active areas

Can detect correlation




# LFP CAPABILITIES

## WHAT IS LFP

- Electrodes can measure cellular spiking activity whose spectrum occur near 1000 Hz as well as the **collective sum of potentials from multiple cells** near the electrode, the local field potential (LFP), usually occurring from 0.1 to 70 Hz
- The Local Field Potential (LFP) is the electric potential recorded in the extracellular space in brain tissue, typically using micro-electrodes (metal, silicon or glass micropipettes). LFPs are recorded in depth, from within the cortical tissue (or other deep brain structures).


# LFP CAPABILITIES

Recording LFP from AO electrodes



# LFP CAPABILITIES

Frequency response



**MACRO TIP**

- Record LFP (0.075Hz-300Hz)
- Macro Stimulation up to 7mA


Zoom In

**MICRO TIP**

- Record Spikes (300Hz-9000Hz)
- Record RAW (0.075Hz-9000Hz)
- Record LFP (0.075Hz-300Hz)
- Micro Stimulation up to 100uA


# LFP CAPABILITIES

Recording LFP from AO electrodes





# LFP CAPABILITIES

RAW data



# LFP CAPABILITIES



Prediction of STN-DBS Electrode Implantation Track in Parkinson's Disease by Using LFP  
*Ilknur Telkes , Joohi Jimenez-Shahed, Ashwin Viswanathan, Aviva Abosch and Nuri F.Ince*

# LFP CAPABILITIES

- Reliable correlations between motor symptoms and components of the LFP power spectrum suggested that LFPs may serve as biomarker for movement disorders
- LFPs are composite signals divided into a number of frequency bands, as follows:
  - 0–3 Hz (delta), 4–7 Hz (theta),
  - 8–12 Hz (alpha), 13–30 Hz (beta), 31–200 Hz (gamma), and >200 Hz (high frequency).

**Table 1.** LFP frequency-based classification system

Frequency band	Range, Hz
Delta	0–3
Theta	4–7
Alpha	8–12
Beta	13–30
Gamma	31–200
High frequency	>200 Hz

# LFP CAPABILITIES

## Neuro Omega system -SETTINGS

- Neuro Omega can record LFP data from Micro and Macro contact of the electrode as well as EEG, EMG and ECOG

### CHANNEL SETTINGS

Channels Settings (Macro)

Name	Chann...	LP F...	HP Fil...	Total ...	S...	SR	Chann...	Contac...	Ref. Co...	Ref. Ap...
Macro...	10005	300....	HW	20.00	On	1375....	LFP	Macro	None	False
Macro...	10389	HW	HW	20.00	On	44000....	RAW	Macro	None	False
Macro...	10006	300....	HW	20.00	On	1375....	LFP	Macro	None	False
Macro...	10390	HW	HW	20.00	On	44000....	RAW	Macro	None	False
Macro...	10007	300....	HW	20.00	On	1375....	LFP	Macro	None	False
Macro...	10391	HW	HW	20.00	On	44000....	RAW	Macro	None	False
Macro...	10008	300....	HW	20.00	On	1375....	LFP	Macro	None	False
Macro...	10392	HW	HW	20.00	On	44000....	RAW	Macro	None	False
Macro...	10009	300....	HW	20.00	Off	1375....	LFP	Macro	None	False
Macro...	10393	HW	HW	20.00	Off	44000....	RAW	Macro	None	False

Change Name   Filter Properties   OK

Default filters are :  
0.07-300Hz

Sampling rate :1375  
Hz

### FILTER PROPERTIES

Filter Properties (Macro)

Macro Contacts:

Ground Contacts

Reference Contacts:

LFP

Turn On

Sampling Rate:

HP (Hz):  Hz

LP (Hz):  Hz

Apply Reference

RAW:

Turn On

Sampling Rate:

Apply Reference

Cancel   Apply   OK

# LFP CAPABILITIES

Neuro Omega system –RELEVANT WINDOWS



Macro LFP

Color Density Spectrogram

