

## Acute Electrodes Order Form

Bill To: (Please provide the billing address)	Deliver To: (Please provide the shipping address)
Name: _____	Name: _____
University: _____	University: _____
Department: _____	Department: _____
Address: _____	Address: _____
Telephone: _____	Telephone: _____
Fax: _____	Fax: _____
Email: _____	Email: _____
Order # _____	_____
VAT # _____	_____

**Refer to next page for a detailed explanation of the different options.**

Number of Boxes <i>10 per Box</i>	Metal	Coating	Shank Diameter $\mu\text{M}$	Total Diameter $\pm 15\% \mu\text{M}$	Length 30-130mm	Taper Angle			Impedance* $\pm 10\% \text{M}\Omega$	Pin <i>Yes or No</i>	Unit Price <i>per Box</i>	Item Total <i>Quantity x Unit Price</i>
						Circle one	choice					
	Tung	Glass	125	250		60	90				153€	
	Tung	Glass	150	310		60	90				192€	
	Tung	Glass	200	390		60	90				192€	
	Tung	Parylene C	125	130	127	S	H	B			135€	
	Tung	Parylene C	250	255	127	S	H	B			135€	
	Plt-Ird	Parylene C	75	80	51	S	H	B			225€	
	Plt-Ird	Parylene C	125	130	76	S	H	B			244€	
	Plt-Ird	Parylene C	250	255	51	S	H	B			225€	
	SS	Parylene C	225	230	76	S	H	B			117€	
	Pur-Ird	Parylene C	100	106	50	S	H	B			315€	
	Pur-Ird	Parylene C + Poly	100	180	50	S	H	B			330€	
	Tung	Parylene C + Poly	125	215	127	S	H	B			183€	
	Tung	Parylene C + Poly	250	355	127	S	H	B			183€	
	Plt-Ird	Parylene C + Poly	75	125	51	S	H	B			259€	
	Plt-Ird	Parylene C + Poly	125	215	76	S	H	B			282€	
	Plt-Ird	Parylene C + Poly	250	355	51	S	H	B			259€	
	SS	Parylene C + Poly	225	355	76	S	H	B			151€	

\* Glass insulated electrodes are available with impedance ranges from 0.5-2.0 M $\Omega$ , except for 390 $\mu\text{M}$  electrodes, 0.5-1.0 M $\Omega$ .

Parylene Coated electrodes are available with impedance ranging from 0.1-2.0 M $\Omega$ .

Subtotal

Shipping based on table to the lower left

**Total not including V.A.T or local taxes**









Shipping and Handling	
1 - 5 boxes	65 €
6 - 10 boxes	85 €
11 - 20 boxes	100 €

**Location and Date** \_\_\_\_\_

**Signature** \_\_\_\_\_

**Please fax this Order to +49 (0)721-2391034**

### Detailed explanation of the choices

Option	Choices	Note
Metal	Tung= Tungsten Plt-Ird= Platinum Iridium SS= Stainless Steel	Some metal choices are not offered with certain insulations.
Coating	Glass	 Glass coated electrodes.
	Parylene C	 Parylene C coated electrodes
	Parylene C with Polyamide	 Polyamide tubing added to enhance stability.
Shank Diameter	Range of choices from 50µM-250µM	Shank diameters are available <b>only</b> as listed in the order form.
Total Diameter	Diameter of glass insulation	Total outer diameter may vary by as much at $\pm 15\% \mu\text{M}$
Length	From 30mm to 130mm depending on coating type	<ul style="list-style-type: none"> <li>- Glass electrodes must be cut to the proper length at the time they are produced. Therefore, customer must specify the length.</li> <li>- Parylene C coated electrodes may be cut to the proper length. The insulation can be burned to expose the wire.</li> </ul>
Taper angel	60 or 90 for glass coating S= Standard H= Heat treated B= Blunt	<div style="display: flex; align-items: center;"> <div style="margin-right: 10px;"> <p>60 or 90</p>  </div> <div style="margin-right: 10px;"> <p>Standard</p>  </div> <div style="margin-right: 10px;"> <p>Heat treated</p>  </div> <div> <p>Blunt</p>  </div> </div>
Pin	0.5 mm standard Male pin	

**Please fax the Order to +49 (0)721-2391034**