



Acute Electrodes Order Form

Bill To: (Please provide the billing address)	Deliver To: (Please provide the shipping address)
Name: _____	Name: _____
University: _____	University: _____
Department: _____	Department: _____
Address: _____	Address: _____
Telephone: _____	Telephone: _____
Fax: _____	Fax: _____
E-mail: _____	E-mail: _____

Refer to next page for a detailed explanation of the different options.

Number of Boxes <i>10 Per Box</i>	Metal	Coating	Shank Diameter μM	Total Diameter $\pm 15\% \mu\text{M}$	Length 30—130mm	Taper Angle Circle one choice	Impedance* $\pm 10\% \text{M}\Omega$	Pin <i>Yes or No</i>	Unit Price <i>Per Box</i>	Item Total <i>Quantity x Unit Price</i>
	Tung	Glass	125	250		60 90			\$200.00	
	Tung	Glass	150	310		60 90			\$250.00	
	Tung	Glass	200	390		60 90			\$250.00	
	Tung	Parylene C	125	130	127	S H B			\$176.00	
	Tung	Parylene C	250	255	127	S H B			\$176.00	
	Plt-Ird	Parylene C	75	80	51	S H B			\$293.00	
	Plt-Ird	Parylene C	125	130	76	S H B			\$318.00	
	Plt-Ird	Parylene C	250	255	51	S H B			\$293.00	
	SS	Parylene C	225	230	76	S H B			\$153.00	
	Pur-Ird	Parylene C	100	106	50	S H B			\$405.00	
	Pur-Ird	Parylene C + Poly	100	180	50	S H B			\$438.00	
	Tung	Parylene C + Poly	125	215	127	S H B			\$239.00	
	Tung	Parylene C + Poly	250	355	127	S H B			\$239.00	
	Plt-Ird	Parylene C + Poly	75	125	51	S H B			\$337.00	
	Plt-Ird	Parylene C + Poly	125	215	76	S H B			\$367.00	
	Plt-Ird	Parylene C + Poly	250	355	51	S H B			\$337.00	
	SS	Parylene C + Poly	225	355	76	S H B			\$197.00	

Shipping and Handling	
1—5 boxes	\$75.00
6—10 boxes	\$95.00
11—20 boxes	\$110.00

Subtotal
Shipping based on table to the left

Total








Credit Card	Purchase Order
Card Type _____	Purchase Order No. _____
Card Number _____ Exp. _____	Purchasing Agent Name _____
Cardholder Name _____	Contact Tel: _____ Fax _____
Billing address _____	Note: - Please fax a copy of the Purchase Order to
City _____ State _____ Zip _____	1(877) 471-2055

Please fax this completed form to (877) 471-2055

* Glass insulated electrodes are available with impedance ranges from 0.5-2.0 M Ω , except for 390 μM electrodes, 0.5-1.0 M Ω .
Parylene Coated electrodes are available with impedance ranging from 0.1-2.0 M Ω .

Acute Electrodes Order Form

Detailed explanation of the choices

Option	Choices	Note
Metal	Tung= Tungsten Plt-Ird= Platinum Iridium SS= Stainless Steel	Some metal choices are not offered with certain insulations.
Coating	Glass	 Glass coated electrodes.
	Parylene C	 Parylene C coated electrodes
	Parylene C with Polyamide	 Polyamide tubing added to enhance stability.
Shank Diameter	Range of choices from 50µM-250µM	Shank diameters are available only as listed in the order form.
Total Diameter	Diameter of glass insulation	Total outer diameter may vary by as much at $\pm 15\% \mu\text{M}$
Length	From 30mm to 130mm depending on coating type	<ul style="list-style-type: none"> - Glass electrodes must be cut to the proper length at the time they are produced. Therefore, customer must specify the length. - Parylene C coated electrodes may be cut to the proper length. The insulation can be burned to expose the wire.
Taper angel	60 or 90 for glass coating S= Standard H= Heat treated B= Blunt	<div style="display: flex; align-items: center;"> <div style="margin-right: 10px;"> <p>60 or 90</p>  </div> <div style="margin-right: 10px;"> <p>Standard</p>  </div> <div style="margin-right: 10px;"> <p>Heat treated</p>  </div> <div> <p>Blunt</p>  </div> </div>
Pin	0.5 mm standard Male pin	