



ALWAYS
INNOVATING
ALWAYS
AHEAD



THE LATEST


Lead Confirm

Alpha Omega has developed a tool to simplify the connectivity between leading DBS leads in the market and Alpha Omega's microelectrode recording systems.

LeadConfirm is used during deep brain stimulation (DBS) surgeries allowing recording and stimulation through the implanted DBS lead.

Key Qualities

- > Connection cables are compatible with leading DBS leads in the market
- > Record from up to 10 contacts of DBS leads with a single Drive Headstage
- > Stimulate through all lead contacts
- > Perform differential recordings between all contacts



Lead
Confirm

